

**SAFETY DATA SHEET**

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 12/13/2006

SDS No : 80-001.5

Date Revised : 06/08/2015

Revision No : 5

**PermaQuick® Bond Resin****SECTION 1 : Identification of the substance/preparation and of the company/undertaking****1.1. Product identifier**

**Product code** : PK/63400  
**Product name** : PermaQuick® Bond Resin  
**Product description** : Dental bonding resin

**1.2. Relevant identified uses of the substance or mixture and uses advised against**

**Relevant identified uses** : Professional Dental Bonding Resin

**1.3. Details of the supplier of the safety data sheet****Manufacturer**

Ultradent Products, Inc.  
 505 W. 10200 S.  
 South Jordan, UT 84095

**Distributor**

Ultradent Products GmbH  
 Am Westhover Berg 30  
 51149 Cologne Germany  
 Email: infoDE@ultradent.com  
**Emergency Phone** : +49(0)2203-35-92-0

**1.4. Emergency telephone number**

**CHEMTREC (NORTH AMERICA)** :(800) 424 - 9300  
**(INTERNATIONAL)** :+1(703) 527 - 3887

**SECTION 2: Hazards identification****2.1. Classification of the substance or mixture****Classification according to Directive 1999/45/EC**

: The mixture has been assessed and/or tested for its physical, health and environmental hazards and the following classification applies.

**Danger symbols** : Xi

**R phrases** : R36/38,R43

**Classification according to Regulation (EC) No 1272/2008 [CLP]**

**Health** : Eye Irritation, Category 2  
 Skin Irritation, Category 2  
 Skin Sensitization, Category 1

**2.2. Label elements****Classification according to Directive 1999/45/EC**

**Hazard pictogram(s)** : Xi



Irritant

**R&S statement(s)** : R36/38: Irritating to eyes and skin.  
 R43: May cause sensitization by skin contact.

**SAFETY DATA SHEET**  
 according to Regulation (EC) No. 453/2010  
 Dental Use

**ULTRADENT**  
 PRODUCTS, INC.

Date Issued : 12/13/2006  
 SDS No : 80-001.5  
 Date Revised : 06/08/2015  
 Revision No : 5

**PermaQuick® Bond Resin**

**Classification according to Regulation (EC) No 1272/2008 [CLP]**

**Hazard pictogram(s)**

:



Exclamation  
mark

**Signal Word**

: WARNING

**Hazard statement(s)**

: H315: Causes skin irritation.  
 H317: May cause an allergic skin reaction.  
 H319: Causes serious eye irritation.

**Precautionary statement(s)**

**Prevention**

: P280: Wear protective gloves/protective clothing/eye protection/face protection.

**Response**

: P305: IF IN EYES:  
 P351: Rinse cautiously with water for several minutes.  
 P315: Get immediate medical advice/attention.  
 P302+P352: IF ON SKIN: Wash with plenty of soap and water.  
 P361: Take off immediately all contaminated clothing.  
 P332+P313: If skin irritation occurs: Get medical advice/attention.  
 P363: Wash contaminated clothing before reuse.

**Storage**

: P273: Avoid release to the environment.

**Disposal**

: P501: Dispose of in compliance with governmental regulation. (EC 1975L0442-20/11/2003)

**2.3. Other hazards**

**Immediate concerns**

: No other adverse human effects.

**SECTION 3: Composition / information on ingredients**

**3.1. Substances**

Not Applicable

**3.2. Mixtures**

**SAFETY DATA SHEET**  
 according to Regulation (EC) No. 453/2010  
 Dental Use



Date Issued : 12/13/2006  
 SDS No : 80-001.5  
 Date Revised : 06/08/2015  
 Revision No : 5

**PermaQuick® Bond Resin**

Chemical Name	CAS	EINECS No.	Wt. %	Classification according to Directive 67/548/EEC	Classification according to Regulation (EC) No 1272/2008 [CLP]
Bis-GMA	1565-94-2	216-367-7	< 20	Xi; R36/37/38; R22	H315; H319; H317
2-hydroxyethyl Methacrylate	868-77-9	212-782-2	< 20	Xi; R43; R36/38	Eye Irr.,Cat. 2; Skin Irr.,Cat. 2; Skin Sens.,Cat. 1; H315; H317; H319
2-dimethylaminoethyl Methacrylate	2867-47-2	220-688-8	< 5	T, Xi; R22; R36/37/38; R43	Acute Tox. (O),Cat. 4; Acute Tox. (D),Cat. 4; Skin Irr.,Cat. 2; Skin Sens.,Cat. 1; Eye Irr.,Cat. 2; H302; H312; H319; H315; H317

For full text of H-statements and R-phrases: see SECTION 16.

**SECTION 4: First aid measures**

**4.1. Description of first aid measures**

- Following eyes** : Immediately flush eyes with plenty of water for at least 15 minutes. Get immediate medical attention.
- Following skin** : Immediately wash skin with soap and plenty of water. Remove contaminated clothing. Get medical attention if symptoms occur. Wash clothing before reuse.
- Following ingestion** : If swallowed in small quantity the product is harmless, otherwise consult a doctor.
- Following inhalation** : No specific treatment is necessary since this material is not likely to be hazardous by inhalation. If exposed to excessive levels of dusts or fumes, remove to fresh air and get medical attention if cough or other symptoms develop.

**4.2. Most important symptoms and effects, both acute and delayed**

**SAFETY DATA SHEET**  
 according to Regulation (EC) No. 453/2010  
 Dental Use



**Date Issued :** 12/13/2006  
**SDS No :** 80-001.5  
**Date Revised :** 06/08/2015  
**Revision No :** 5

**PermaQuick® Bond Resin**

**Eyes** : Causes eye irritation.  
**Skin** : May cause skin irritation and sensitization.  
**Ingestion** : May be harmful if swallowed in large quantities.  
**Inhalation** : None expected for this product.

**4.3. Indication of any immediate medical attention and special treatment needed**

**Notes to physician** : Not Applicable

**SECTION 5: Fire fighting measures**

**5.1. Extinguishing media**

**Extinguishing media** : Please see Fire Fighting Equipment under Section 5.3.

**5.2. Special hazards arising from the substance or mixture**

**Hazardous combustion products** : None known.  
**Explosion hazards** : Not Determined  
**Fire explosion** : Not Determined  
**Sensitive to static discharge** : Not Determined  
**Sensitivity to impact** : Not Determined

**5.3. Advice for firefighters**

**Fire fighting procedures** : General: Evacuate all personnel; use protective equipment for fire-fighting. Use self-contained breathing apparatus when the product is involved in fire.  
**Fire fighting equipment** : Foam, dry chemical, carbon dioxide (CO<sub>2</sub>).

**SECTION 6: Accidental release measures**

**6.1. Personal precautions, protective equipment and emergency procedures**

**General procedures** : Refer to Section 8 for Personal Protective Equipment.

**6.2. Environmental precautions**

**Water spill** : Do not allow to enter sewers or drains that may lead to waterways.

**6.3. Methods and material for containment and cleaning up**

**Small spill** : Pick up wash liquid with additional absorbent and place in a disposable container.  
**Large spill** : Collect in tightly closed containers.

**6.4. Reference to other sections**

**Reference to other sections** : Not Applicable

**SECTION 7: Handling and storage**

**7.1. Precautions for safe handling**

**Handling** : Use suitable protective equipment.  
**Storage** : See product labeling.

**SAFETY DATA SHEET**

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 12/13/2006

SDS No : 80-001.5

Date Revised : 06/08/2015

Revision No : 5

**PermaQuick® Bond Resin****7.2. Conditions for safe storage, including any incompatibilities**

Shelf life : See product labeling

**7.3. Specific end use(s)**

Specific end use(s) : Professional Dental Bonding Resin

**SECTION 8: Exposure controls / personal protection****8.1. Control parameters**

Control parameters : Not Determined

**8.2. Exposure controls**

Eye/face protection : Wear eye protection

Skin protection : Wear suitable protective clothing and gloves.

Respiratory protection : Good general ventilation should be sufficient to control airborne levels.

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state : Viscous substance

Colour : Beige

Odour : Acrylic

**9.2. Other information**

Percent volatile : Not Determined

**SECTION 10: Stability and reactivity****10.1. Reactivity**

Reactivity : Stable

**10.2. Chemical stability**

Chemical stability : Stable when stored and handled under recommended conditions.

**10.3. Possibility of hazardous reactions**

Hazardous polymerization : Stable when stored and handled under recommended conditions.

**10.4. Conditions to avoid**

Conditions to avoid : Heat, flames, ignition sources, direct light, and incompatibles.

**10.5. Incompatible materials**

Incompatible materials : Not Determined

**10.6. Hazardous decomposition products**

Hazardous decomposition products : None known.

**SECTION 11: Toxicological information**

**SAFETY DATA SHEET**

according to Regulation (EC) No. 453/2010

Dental Use



**Date Issued :** 12/13/2006  
**SDS No :** 80-001.5  
**Date Revised :** 06/08/2015  
**Revision No :** 5

**PermaQuick® Bond Resin****11.1. Information on toxicological effects****Acute**

**Notes** : Device is biocompatible when used as directed by dental professionals per ISO 10993-1.

**SECTION 12: Ecological information****12.1. Toxicity****Aquatic toxicity (acute)**

**96-hour LC<sub>50</sub>** : Not Determined  
**48-hour EC<sub>50</sub>** : Not Determined  
**96-hour EC<sub>50</sub>** : Not Determined

**12.2. Persistence and degradability**

**Persistence and degradability** : Not Determined

**12.3. Bioaccumulative potential**

**Bioaccumulative potential** : Not Determined

**12.4. Mobility in soil**

**Mobility in soil** : Not Determined

**12.5. Results of PBT and vPvB assessment**

**Results of PBT and vPvB assessment** : Not Determined

**12.6. Other adverse effects**

**Environmental data** : Do not flush to sewer.

**SECTION 13: Disposal considerations****13.1. Waste treatment methods**

**Disposal method** : Dispose of in compliance with governmental regulation.(EC 1975L0442-20/11/2003)

**SECTION 14: Transport information****14.1. UN number**

**UN number** : N/A

**14.2. UN proper shipping name**

**UN proper shipping name** : Not Dangerous Goods

**14.3. Transport hazard class(es)**

**Hazard classification** : N/A

**14.4. Packing group**

**Packing group** : N/A

**14.5. Environmental hazards**

**SAFETY DATA SHEET**

according to Regulation (EC) No. 453/2010

Dental Use



Date Issued : 12/13/2006

SDS No : 80-001.5

Date Revised : 06/08/2015

Revision No : 5

**PermaQuick® Bond Resin**

Marine pollutant #1 : N/A

**14.6. Special precautions for user**

ADR - road : Not Dangerous Goods

RID - rail : Not Dangerous Goods

IMDG - sea : Not Dangerous Goods

IATA - air : Not Dangerous Goods

**14.7. Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code**

Transport in bulk : Not Dangerous Goods

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture**

RoHS : Please refer to Medical Devices Directive 93/42/EEC

**15.2. Chemical safety assessment**

Chemical safety assessment : See Section 11

Additional information : This product is not considered dangerous for transportation

**SECTION 16: Other information**

**Relevant R-phrases and/or H-statements (number and full text)** : R22: Harmful if swallowed.  
 R36/37/38: Irritating to eyes, respiratory system and skin.  
 R36/38: Irritating to eyes and skin.  
 R43: May cause sensitization by skin contact.  
 Acute Tox. (D), Cat. 4: Acute Toxicity (Dermal), Category 4  
 Acute Tox. (O), Cat. 4: Acute Toxicity (Oral), Category 4  
 Eye Irr., Cat. 2: Eye Irritation, Category 2  
 Skin Irr., Cat. 2: Skin Irritation, Category 2  
 Skin Sens., Cat. 1: Skin Sensitization, Category 1  
 H302: Harmful if swallowed.  
 H312: Harmful in contact with skin.  
 H315: Causes skin irritation.  
 H317: May cause an allergic skin reaction.  
 H319: Causes serious eye irritation.

Prepared by : Anu Kattoju

Revision summary : This SDS replaces the 04/06/2015 SDS. Revised: **Section 14:**  
**Section 15:** COMMENTS.

General statements : N/A= Not Applicable

**Manufacturer disclaimer** : FOR DENTAL USE ONLY: Use as directed. The information and recommendations are taken from sources (raw material SDS(s) and manufacturer's knowledge) believed to be accurate; however, the manufacturer, makes no warranty with respect to the accuracy of the information or the suitability of the recommendation and assumes no liability to any user thereof. Each user should review these recommendations in the specific context of the intended use and

# SAFETY DATA SHEET

according to Regulation (EC) No. 453/2010

Dental Use



**Date Issued :** 12/13/2006

**SDS No :** 80-001.5

**Date Revised :** 06/08/2015

**Revision No :** 5

## PermaQuick® Bond Resin

---

determine whether they are appropriate.